



## The use of HiSPECT in the investigation of the 5-HT<sub>2A</sub> receptor and DAT status in the Mongolian gerbil.

### 01 | Background

Stereotypies are ethologically defined as repetitive, invariant behaviours, apparently without function (Ödberg 1978). A number of studies show that stereotypies in animals are caused by a chronic frustration resulting from a stimulus-poor captive environment (for a review see Mason 1991). Humans also display stereotypies, in relation to mental conditions such as schizophrenia, bipolar disorder, dementia (Mendez et al 2005), autism (Lewis & Bodfish 1998),... or during moments of inner tension (Soussignan et al 1988). The involvement of the dopaminergic and serotonergic system in the occurrence of drug-induced and even spontaneous stereotypies has been widely accepted (reviewed in Cabib 2006, Pogorelov et al. 2005, Meers & Ödberg 2005).

Therefore, gaining more insight in the serotonergic and dopaminergic status in animal models of stereotypy is useful for the future study of animal stereotypies, but can also point out which markers are additionally useful in the study of certain human psychiatric conditions.

### 02 | Aims

- 1) Investigate the feasibility of using multi-pinhole high resolution SPECT imaging in Mongolian gerbils and verify whether a conventional perfusion tracer will allow us to anatomically identify gerbil brain structures.
- 2) Investigate the appropriateness of the radioligand for DAT imaging.
- 3) Investigate the appropriateness of the radioligand for 5-HT<sub>2A</sub> imaging.

### 03 | Materials and Methods

#### Animals

- 16 group-housed (2-3) Mongolian gerbils (355 ± 65 days)
- Makrolon IV cages, standard laboratory conditions
- Reversed 14:10 light-dark cycle (light: 19:00 - 9:00)
- Environmental Enrichment: chewing blocks, paper tissue and hay
- General anaesthesia during scans: (4% isoflurane for induction, 1,8% for maintenance)

#### High-resolution SPECT

- HiSPECT (Bioscan) using conventional triple head gamma camera with multi-pinhole collimation (6 holes with 1,5mm diameter)
- Step & Shoot: 36°, 10 steps, 120 s per step, 256 x 256 matrix
- Image reconstruction using a dedicated OSEM algorithm (Scivis)
- IV injection of tracer into femoral vein



#### References:

Cabib S 2006. The neurobiology of stereotypy II: the role of stress. In: *Stereotypic Animal Behaviour*, eds. Mason G & J Rushen, p227-255. Lewis MH & Bodfish JW 1998. Repetitive behavior disorders in autism. *Mental Retardation and Developmental Disabilities Research Reviews* 4(2):80-89. Ödberg FO 1978. Abnormal behaviours (Stereotypies). In: *Proceedings of the 1st World Congress on Ethology Applied in Zootechnics* (Editorial Garsi), p. 475-480. Madrid: Industrias Graficas España. Mason GJ 1991. Stereotypies: a critical review. *Animal Behaviour* 41(6): 1015-1037. Meers L & FO Ödberg 2005. Paradoxical rate-dependent effect of fluoxetine on captivity-induced stereotypies in bank voles. *Progress in Neuro-Psychopharmacology* 29:962-969. Mendez MF, Shapira JS & JL Miller 2005. Stereotypical movements and frontotemporal dementia. *Movement disorders* 20(6): 742-745. Pogorelov VM, Rodriguez RM, Insko ML, Caron MG & WC Wetsel 2005. Novelty seeking and stereotypic activation of behavior in mice with disruption of DAT1 gene. *Neuropsychopharmacology* 30: 1818-1831. Soussignan R, Koch P & H Montagner 1988. Relationship between the general motor-activity of children at school and the duration of school tasks which are imposed on them. *Comptes rendus de l'Academie des Sciences serie III - sciences de la vie* 306:139-142.

#### Ghent University, Belgium,

- Moons CPH Department of Nutrition, Genetics & Ethology, Faculty of Veterinary Medicine
- Peremans K Medical imaging of domestic animals, Faculty of Veterinary Medicine
- Audenaert K Department of Psychiatry and Medical Psychology, Faculty of Medicine and Health Sciences
- Dobbeleir A Medical imaging of domestic animals, Faculty of Veterinary Medicine
- Vandermeulen E Medical imaging of domestic animals, Faculty of Veterinary Medicine
- Vermeire S Medical imaging of domestic animals, Faculty of Veterinary Medicine
- De Vos F Department of Pharmaceutical Analysis, Faculty of Pharmaceutical Sciences
- Soenen E Internal Medicine and Clinical Biology of Large Animals, Faculty of Veterinary Medicine
- Van Bree H Medical imaging of domestic animals, Faculty of Veterinary Medicine
- Ödberg FO Department of Nutrition, Genetics & Ethology, Faculty of Veterinary Medicine


### 04 | Results

#### 1) Perfusion data for anatomical identification

- 4 gerbils
- 24,42 ± 5,92 MBq <sup>99</sup>Tc-ethylcysteinate Dimer
- Scanning 15 min after injection
- Perfusion image shows that spinal cord, cerebellum and cortical structures are easily distinguished.

#### 2) DAT

- 5 gerbils (1 for scanning optimum confirmation)
- 44,33 ± 11,66 MBq <sup>123</sup>I-FP CIT
- Scanning optimum: 4 ½ hours after injection
- Specific binding to the DAT sites observed by the regional distribution of the radioactivity in the brain which correlated well with the known presence of DAT sites (basal ganglia, white arrows).




#### 3) 5-HT<sub>2A</sub>

- 4 gerbils
- 22,2 ± 3,7 MBq <sup>123</sup>I-R91-150
- 90 min after injection

#### Two problems were encountered:

- First, the images were initially of poor quality. This problem was solved using Cyclosporine A (50 mg/kg) to block the P-glycoprotein pump in the blood-brain barrier.



- Second, when determining the scanning optimum, the cortex/cerebellum ratio of activity reached a rather unstable plateau. Further studies are being conducted looking into the second problem.

### 05 | Conclusion

- 1) It is feasible to use the conventional ECD tracer to perform perfusion imaging for anatomical identification.
- 2) The FP CIT tracer is appropriate to monitor the DAT status in the gerbil.
- 3) The results for the 5-HT<sub>2A</sub> receptor were less straightforward, but nonetheless promising and are being investigated further.